

Date	Days Post-Transduction	Days Post-Treatment	Selection Agent	Drug	Concentration (nM)	Actual concentration replated (nM)	Replicate	Plate # /Format	Cell Count (millions)	Cell Count Replated	# gDNA cell pellets	gDNA pellet cell count (millions)	RFP+ %	Notes
9/30/24	0		Puromycin (5 ug/mL)	NONE	N/A	N/A	1	10x 6-well plates	120	N/A	N/A	N/A	38.40%	
9/30/24	0		Puromycin (5 ug/mL)	NONE	N/A	N/A	2	10x 6-well plates	120	N/A	N/A	N/A	35.70%	
9/30/24	0		Puromycin (5 ug/mL)	NONE	N/A	N/A	3	10x 6-well plates	120	N/A	N/A	N/A	37.70%	
10/1/24	1		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	1	15x 15-cm plates	N/A	N/A	N/A	N/A	N/A	
10/1/24	1		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	2	15x 15-cm plates	N/A	N/A	N/A	N/A	N/A	
10/1/24	1		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	3	15x 15-cm plates	N/A	N/A	N/A	N/A	N/A	
10/4/24	4		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	1	15x 15-cm plates	189	70	1	120	N/A	
10/4/24	4		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	2	15x 15-cm plates	187	70	1	120	N/A	
10/4/24	4		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	3	15x 15-cm plates	186.5	70	1	120	N/A	
10/7/24	7		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	1	13x 15-cm plates	118	36	4	3	79.80%	Additional cell pellet taken with remaining cells = ~40-50 MM
10/7/24	7		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	2	13x 15-cm plates	124	36	4	3	81.90%	Additional cell pellet taken with remaining cells = ~40-50 MM
10/7/24	7		Puromycin (5 ug/mL) Blasticidin (10 ug/mL)	NONE	N/A	N/A	3	13x 15-cm plates	114	36	4	3	80.50%	Additional cell pellet taken with remaining cells = ~40-50 MM
10/8/24	8	0	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	Cell counts are based on cell counts at time of plating 24 hours prior. Drug is added to cells at this time-point, after cells adhere.
10/8/24	8	0	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32	1000	1000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32	1000	1000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32	1000	1000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32N	5000	5000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32N	5000	5000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KI-CDK9d-32N	5000	5000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KB-0742	1500	1500	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KB-0742	1500	1500	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KB-0742	1500	1500	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KB-0742	5000	5000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KB-0742	5000	5000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	KB-0742	5000	5000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SY-5609	1000	1000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SY-5609	1000	1000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SY-5609	1000	1000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SY-5609	2500	2500	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SY-5609	2500	2500	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SY-5609	2500	2500	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	Senexin B	2000	2000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	Senexin B	2000	2000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	Senexin B	2000	2000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	Senexin B	15000	15000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	Senexin B	15000	15000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	Senexin B	15000	15000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SEL120	4000	4000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SEL120	4000	4000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SEL120	4000	4000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SEL120	8000	8000	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SEL120	8000	8000	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	SEL120	8000	8000	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	DMSO	1%	1%	1	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	DMSO	1%	1%	2	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/8/24	8	0	NONE	DMSO	1%	1%	3	1x 15-cm plate	2.75	N/A	N/A	N/A	N/A	
10/11/24	11	3	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	6.93	2.75	1	3	N/A	
10/11/24	11	3	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	5.94	2.75	1	3	N/A	
10/11/24	11	3	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	6.64	2.75	1	3	N/A	
10/11/24	11	3	NONE	KI-CDK9d-32	1000	1000	1	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KI-CDK9d-32	1000	1000	2	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KI-CDK9d-32	1000	1000	3	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	7.57	2.75	1	3	N/A	
10/11/24	11	3	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	6.72	2.75	1	3	N/A	
10/11/24	11	3	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	7.94	2.75	1	3	N/A	
10/11/24	11	3	NONE	KI-CDK9d-32N	5000	5000	1	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KI-CDK9d-32N	5000	5000	2	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KI-CDK9d-32N	5000	5000	3	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KB-0742	1500	1500	1	1x 15-cm plate	6.98	2.75	1	3	N/A	
10/11/24	11	3	NONE	KB-0742	1500	1500	2	1x 15-cm plate	5.71	2.75	1	3	N/A	
10/11/24	11	3	NONE	KB-0742	1500	1500	3	1x 15-cm plate	6.46	2.75	1	3	N/A	
10/11/24	11	3	NONE	KB-0742	5000	5000	1	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KB-0742	5000	5000	2	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	KB-0742	5000	5000	3	1x 15-cm plate	N/A				N/A	Media change with fresh drug
10/11/24	11	3	NONE	SY-5609	1000	1000	1	1x 15-cm plate	1.8	2.75	N/A	N/A	N/A	No pellet taken; replated all cells
10/11/24	11	3	NONE	SY-5609	1000	1000	2	1x 15-cm plate	2.43	2.75	N/A	N/A	N/A	No pellet taken; replated all cells
10/11/24	11	3	NONE	SY-5609	1000	1000	3	1x 15-cm plate	1.61	2.75	N/A	N/A	N/A	No pellet taken; replated all cells
10/11/24	11	3	NONE	SY-5609	2500	2500	1	1x 15-cm plate	N/A				N/A	
10/11/24	11	3	NONE	SY-5609	2500	2500	2	1x 15-cm plate	N/A				N/A	
10/11/24	11	3	NONE	SY-5609	2500	2500	3	1x 15-cm plate	N/A				N/A	
10/11/24	11	3	NONE	Senexin B	2000	2000	1	1x 15-cm plate	8.72	2.75	1	3	N/A	
10/11/24	11	3	NONE	Senexin B	2000	2000	2	1x 15-cm plate	8.62	2.75	1	3	N/A	
10/11/24	11	3	NONE	Senexin B	2000	2000	3	1x 15-cm plate	7.26	2.75	1	3	N/A	

10/11/24	11	3	NONE	Senexin B	15000	15000	1	1x 15-cm plate	N/A					N/A	Media change with fresh drug
10/11/24	11	3	NONE	Senexin B	15000	15000	2	1x 15-cm plate	N/A					N/A	Media change with fresh drug
10/11/24	11	3	NONE	Senexin B	15000	15000	3	1x 15-cm plate	N/A					N/A	Media change with fresh drug
10/11/24	11	3	NONE	SEL120	4000	4000	1	1x 15-cm plate	8.35	2.75	1		3	N/A	
10/11/24	11	3	NONE	SEL120	4000	4000	2	1x 15-cm plate	8.62	2.75	1		3	N/A	
10/11/24	11	3	NONE	SEL120	4000	4000	3	1x 15-cm plate	7.26	2.75	1		3	N/A	
10/11/24	11	3	NONE	SEL120	8000	8000	1	1x 15-cm plate	7.01	2.75	1		3	N/A	
10/11/24	11	3	NONE	SEL120	8000	8000	2	1x 15-cm plate	7.33	2.75	1		3	N/A	
10/11/24	11	3	NONE	SEL120	8000	8000	3	1x 15-cm plate	6.92	2.75	1		3	N/A	
10/11/24	11	3	NONE	DMSO	1%	1%	1	1x 15-cm plate	10.4	2.75	2		3	N/A	
10/11/24	11	3	NONE	DMSO	1%	1%	2	1x 15-cm plate	9	2.75	2		3	N/A	
10/11/24	11	3	NONE	DMSO	1%	1%	3	1x 15-cm plate	9.86	2.75	2		3	N/A	
10/14/24	14	6	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	5.34	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	5.31	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	5.92	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KI-CDK9d-32	1000	0	1	1x 15-cm plate	1.78	1.78				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KI-CDK9d-32	1000	0	2	1x 15-cm plate	1.84	1.84				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KI-CDK9d-32	1000	0	3	1x 15-cm plate	1.9	1.9				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	6.14	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	5.01	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	6.23	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KI-CDK9d-32N	5000	0	1	1x 15-cm plate	1.08	1.08				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KI-CDK9d-32N	5000	0	2	1x 15-cm plate	1.19	1.19				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KI-CDK9d-32N	5000	0	3	1x 15-cm plate	1.29	1.29				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KB-0742	1500	1500	1	1x 15-cm plate	5.03	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KB-0742	1500	1500	2	1x 15-cm plate	5.56	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KB-0742	1500	1500	3	1x 15-cm plate	5.39	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	KB-0742	5000	0	1	1x 15-cm plate	1.51	1.51				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KB-0742	5000	0	2	1x 15-cm plate	1.22	1.22				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	KB-0742	5000	0	3	1x 15-cm plate	1.34	1.34				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	SY-5609	1000	0	1	1x 15-cm plate	1.32	1.32				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	SY-5609	1000	0	2	1x 15-cm plate	1.59	1.59				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	SY-5609	1000	0	3	1x 15-cm plate	1.48	1.48				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	SY-5609	2500	0	1	1x 15-cm plate	1.59	1.59				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	SY-5609	2500	0	2	1x 15-cm plate	1.28	1.28				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	SY-5609	2500	0	3	1x 15-cm plate	1.16	1.16				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	Senexin B	2000	2000	1	1x 15-cm plate	4.61	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	Senexin B	2000	2000	2	1x 15-cm plate	4.66	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	Senexin B	2000	2000	3	1x 15-cm plate	5.48	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	Senexin B	15000	0	1	1x 15-cm plate	2.41	2.41				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	Senexin B	15000	0	2	1x 15-cm plate	2.49	2.49				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	Senexin B	15000	0	3	1x 15-cm plate	2.63	2.63				N/A	Initially plated cells in drug on 10/14. Replace with media with no drug on 10/15 morning.
10/14/24	14	6	NONE	SEL120	4000	4000	1	1x 15-cm plate	5.81	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	SEL120	4000	4000	2	1x 15-cm plate	5.46	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	SEL120	4000	4000	3	1x 15-cm plate	6.25	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	SEL120	8000	8000	1	1x 15-cm plate	5.27	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	SEL120	8000	8000	2	1x 15-cm plate	4.38	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	SEL120	8000	8000	3	1x 15-cm plate	4.77	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	DMSO	1%	1%	1	1x 15-cm plate	7.48	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	DMSO	1%	1%	2	1x 15-cm plate	6.66	2.75	1		2.75	N/A	
10/14/24	14	6	NONE	DMSO	1%	1%	3	1x 15-cm plate	5.55	2.75	1		2.75	N/A	
10/17/24	17	9	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	5.67	2.75	1		2.75	N/A	
10/17/24	17	9	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	5.49	2.75	1		2.75	N/A	
10/17/24	17	9	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	5.95	2.75	1		2.75	N/A	
10/17/24	17	9	NONE	KI-CDK9d-32	1000	0	1	1x 15-cm plate	N/A					N/A	
10/17/24	17	9	NONE	KI-CDK9d-32	1000	0	2	1x 15-cm plate	N/A					N/A	
10/17/24	17	9	NONE	KI-CDK9d-32	1000	0	3	1x 15-cm plate	N/A					N/A	
10/17/24	17	9	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	6.14	2.75	1		2.75	N/A	
10/17/24	17	9	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	6.19	2.75	1		2.75	N/A	

10/17/24	17	9	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	5.9	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	KI-CDK9d-32N	5000	0	1	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	KI-CDK9d-32N	5000	0	2	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	KI-CDK9d-32N	5000	0	3	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	KB-0742	1500	1500	1	1x 15-cm plate	6.3	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	KB-0742	1500	1500	2	1x 15-cm plate	5.3	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	KB-0742	1500	1500	3	1x 15-cm plate	6.68	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	KB-0742	5000	0	1	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	KB-0742	5000	0	2	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	KB-0742	5000	0	3	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	SY-5609	1000	0	1	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	SY-5609	1000	0	2	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	SY-5609	1000	0	3	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	SY-5609	2500	0	1	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	SY-5609	2500	0	2	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	SY-5609	2500	0	3	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	Senexin B	2000	2000	1	1x 15-cm plate	6.71	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	Senexin B	2000	2000	2	1x 15-cm plate	5.6	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	Senexin B	2000	2000	3	1x 15-cm plate	6.05	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	Senexin B	15000	0	1	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	Senexin B	15000	0	2	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	Senexin B	15000	0	3	1x 15-cm plate	N/A				N/A	
10/17/24	17	9	NONE	SEL120	4000	4000	1	1x 15-cm plate	7.17	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	SEL120	4000	4000	2	1x 15-cm plate	7.02	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	SEL120	4000	4000	3	1x 15-cm plate	6.47	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	SEL120	8000	8000	1	1x 15-cm plate	2.8	2.8			N/A	
10/17/24	17	9	NONE	SEL120	8000	8000	2	1x 15-cm plate	2.18	2.18			N/A	
10/17/24	17	9	NONE	SEL120	8000	8000	3	1x 15-cm plate	3.49	3.49			N/A	
10/17/24	17	9	NONE	DMSO	1%	1%	1	1x 15-cm plate	7.17	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	DMSO	1%	1%	2	1x 15-cm plate	7.17	2.75	1	2.75	N/A	
10/17/24	17	9	NONE	DMSO	1%	1%	3	1x 15-cm plate	6.98	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	4.85	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	4.69	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	4.49	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KI-CDK9d-32	1000	0	1	1x 15-cm plate	3.66	3.66			N/A	
10/20/24	20	12	NONE	KI-CDK9d-32	1000	0	2	1x 15-cm plate	3.15	3.15			N/A	
10/20/24	20	12	NONE	KI-CDK9d-32	1000	0	3	1x 15-cm plate	3.29	3.29			N/A	
10/20/24	20	12	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	5.07	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	5.81	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	5.99	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KI-CDK9d-32N	5000	0	1	1x 15-cm plate	3.36	3.36			N/A	
10/20/24	20	12	NONE	KI-CDK9d-32N	5000	0	2	1x 15-cm plate	2.29	2.29			N/A	
10/20/24	20	12	NONE	KI-CDK9d-32N	5000	0	3	1x 15-cm plate	3.28	3.28			N/A	
10/20/24	20	12	NONE	KB-0742	1500	1500	1	1x 15-cm plate	4.77	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KB-0742	1500	1500	2	1x 15-cm plate	4.95	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KB-0742	1500	1500	3	1x 15-cm plate	4.07	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	KB-0742	5000	0	1	1x 15-cm plate	0.115	0.115			N/A	
10/20/24	20	12	NONE	KB-0742	5000	0	2	1x 15-cm plate	0.183	0.183			N/A	
10/20/24	20	12	NONE	KB-0742	5000	0	3	1x 15-cm plate	0.0994	0.0994			N/A	
10/20/24	20	12	NONE	SY-5609	1000	0	1	1x 15-cm plate	0.816	0.816			N/A	
10/20/24	20	12	NONE	SY-5609	1000	0	2	1x 15-cm plate	0.821	0.821			N/A	
10/20/24	20	12	NONE	SY-5609	1000	0	3	1x 15-cm plate	0.649	0.649			N/A	
10/20/24	20	12	NONE	SY-5609	2500	0	1	1x 15-cm plate	0.685	0.685			N/A	
10/20/24	20	12	NONE	SY-5609	2500	0	2	1x 15-cm plate	0.659	0.659			N/A	
10/20/24	20	12	NONE	SY-5609	2500	0	3	1x 15-cm plate	0.602	0.602			N/A	
10/20/24	20	12	NONE	Senexin B	2000	2000	1	1x 15-cm plate	4.89	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	Senexin B	2000	2000	2	1x 15-cm plate	4.92	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	Senexin B	2000	2000	3	1x 15-cm plate	5.22	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	Senexin B	15000	4000	1	1x 15-cm plate	5.07	2.75	1	2.75	N/A	Replated in new concentration (4000 nM)
10/20/24	20	12	NONE	Senexin B	15000	4000	2	1x 15-cm plate	5.49	2.75	1	2.75	N/A	Replated in new concentration (4000 nM)
10/20/24	20	12	NONE	Senexin B	15000	4000	3	1x 15-cm plate	5.19	2.75	1	2.75	N/A	Replated in new concentration (4000 nM)
10/20/24	20	12	NONE	SEL120	4000	4000	1	1x 15-cm plate	6.72	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	SEL120	4000	4000	2	1x 15-cm plate	6.83	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	SEL120	4000	4000	3	1x 15-cm plate	6.31	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	SEL120	8000	8000	1	1x 15-cm plate	4.55	4.55			N/A	
10/20/24	20	12	NONE	SEL120	8000	8000	2	1x 15-cm plate	1.25	1.25			N/A	
10/20/24	20	12	NONE	SEL120	8000	8000	3	1x 15-cm plate	3.46	3.46			N/A	
10/20/24	20	12	NONE	DMSO	1%	1%	1	1x 15-cm plate	6.73	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	DMSO	1%	1%	2	1x 15-cm plate	7.06	2.75	1	2.75	N/A	
10/20/24	20	12	NONE	DMSO	1%	1%	3	1x 15-cm plate	6.26	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	5.14	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	5.05	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	5.44	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32	1000	200	1	1x 15-cm plate	7.09	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32	1000	200	2	1x 15-cm plate	7.17	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32	1000	200	3	1x 15-cm plate	7.55	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	4.99	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	5.92	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	5.52	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32N	5000	2000	1	1x 15-cm plate	6.51	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32N	5000	2000	2	1x 15-cm plate	9.02	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KI-CDK9d-32N	5000	2000	3	1x 15-cm plate	8.74	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KB-0742	1500	1500	1	1x 15-cm plate	4.79	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KB-0742	1500	1500	2	1x 15-cm plate	4.88	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KB-0742	1500	1500	3	1x 15-cm plate	4.66	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	KB-0742	5000	0	1	1x 15-cm plate						Terminated these cells. No growth
10/23/24	23	15	NONE	KB-0742	5000	0	2	1x 15-cm plate						Terminated these cells. No growth
10/23/24	23	15	NONE	KB-0742	5000	0	3	1x 15-cm plate						Terminated these cells. No growth
10/23/24	23	15	NONE	SY-5609	1000	0	1	1x 15-cm plate	0.356	0.356			N/A	
10/23/24	23	15	NONE	SY-5609	1000	0	2	1x 15-cm plate	0.513	0.513			N/A	
10/23/24	23	15	NONE	SY-5609	1000	0	3	1x 15-cm plate	0.309	0.309			N/A	
10/23/24	23	15	NONE	SY-5609	2500	0	1	1x 15-cm plate	0.403	0.403			N/A	
10/23/24	23	15	NONE	SY-5609	2500	0	2	1x 15-cm plate	0.251	0.251			N/A	
10/23/24	23	15	NONE	SY-5609	2500	0	3	1x 15-cm plate	0.398	0.398			N/A	
10/23/24	23	15	NONE	Senexin B	2000	2000	1	1x 15-cm plate	3.67	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	Senexin B	2000	2000	2	1x 15-cm plate	3.81	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	Senexin B	2000	2000	3	1x 15-cm plate	5.34	2.75	1	2.75	N/A	
10/23/24	23	15	NONE	Senexin B	15000	0	1	1x 15-cm plate	3.36	3.36			N/A	

10/23/24	23	15	NONE	Senexin B	15000	0	2	1x 15-cm plate	3.36	3.36				N/A	
10/23/24	23	15	NONE	Senexin B	15000	0	3	1x 15-cm plate	3.15	3.15				N/A	
10/23/24	23	15	NONE	SEL120	4000	4000	1	1x 15-cm plate	5.02	2.75	1		2.75	N/A	
10/23/24	23	15	NONE	SEL120	4000	4000	2	1x 15-cm plate	5.69	2.75	1		2.75	N/A	
10/23/24	23	15	NONE	SEL120	4000	4000	3	1x 15-cm plate	5.11	2.75	1		2.75	N/A	
10/23/24	23	15	NONE	SEL120	8000	0	1	1x 15-cm plate	4.15	4.15				N/A	
10/23/24	23	15	NONE	SEL120	8000	0	2	1x 15-cm plate	0.141	0.141				N/A	
10/23/24	23	15	NONE	SEL120	8000	0	3	1x 15-cm plate	1.2	1.2				N/A	
10/23/24	23	15	NONE	DMSO	1%	1%	1	1x 15-cm plate	6.45	2.75	1		2.75	N/A	
10/23/24	23	15	NONE	DMSO	1%	1%	2	1x 15-cm plate	6.67	2.75	1		2.75	N/A	
10/23/24	23	15	NONE	DMSO	1%	1%	3	1x 15-cm plate	7.75	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	4.95	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	5.04	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	5.29	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32	1000	200	1	1x 15-cm plate	4.3	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32	1000	200	2	1x 15-cm plate	4.41	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32	1000	200	3	1x 15-cm plate	4.23	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	5.67	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	5.19	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	6.04	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32N	5000	2000	1	1x 15-cm plate	4.22	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32N	5000	2000	2	1x 15-cm plate	4.17	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KI-CDK9d-32N	5000	2000	3	1x 15-cm plate	4.05	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KB-0742	1500	1500	1	1x 15-cm plate	4.17	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KB-0742	1500	1500	2	1x 15-cm plate	4.8	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KB-0742	1500	1500	3	1x 15-cm plate	4.36	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	KB-0742	5000	0	1	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	KB-0742	5000	0	2	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	KB-0742	5000	0	3	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	SY-5609	1000	0	1	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	SY-5609	1000	0	2	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	SY-5609	1000	0	3	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	SY-5609	2500	0	1	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	SY-5609	2500	0	2	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	SY-5609	2500	0	3	1x 15-cm plate						N/A	
10/26/24	26	18	NONE	Senexin B	2000	2000	1	1x 15-cm plate	4.68	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	Senexin B	2000	2000	2	1x 15-cm plate	4.8	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	Senexin B	2000	2000	3	1x 15-cm plate	5.52	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	Senexin B	15000	4000	1	1x 15-cm plate	6.33	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	Senexin B	15000	4000	2	1x 15-cm plate	6.15	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	Senexin B	15000	4000	3	1x 15-cm plate	6.16	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	SEL120	4000	4000	1	1x 15-cm plate	4.62	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	SEL120	4000	4000	2	1x 15-cm plate	4.96	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	SEL120	4000	4000	3	1x 15-cm plate	3.66	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	SEL120	8000	0	1	1x 15-cm plate	7.13	2.75				N/A	
10/26/24	26	18	NONE	SEL120	8000	0	2	1x 15-cm plate	0.00471					N/A	Terminated this replicate. NO cells
10/26/24	26	18	NONE	SEL120	8000	0	3	1x 15-cm plate	0.978	0.978				N/A	
10/26/24	26	18	NONE	DMSO	1%	1%	1	1x 15-cm plate	7.87	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	DMSO	1%	1%	2	1x 15-cm plate	7.1	2.75	1		2.75	N/A	
10/26/24	26	18	NONE	DMSO	1%	1%	3	1x 15-cm plate	7.49	2.75	1		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32	100	100	1	1x 15-cm plate	5.57		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32	100	100	2	1x 15-cm plate	5.86		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32	100	100	3	1x 15-cm plate	6.62		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32	1000	200	1	1x 15-cm plate	6.73		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32	1000	200	2	1x 15-cm plate	6.21		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32	1000	200	3	1x 15-cm plate	5.36		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32N	1250	1250	1	1x 15-cm plate	4.84		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32N	1250	1250	2	1x 15-cm plate	5.72		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32N	1250	1250	3	1x 15-cm plate	5.34		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32N	5000	2000	1	1x 15-cm plate	6.69		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32N	5000	2000	2	1x 15-cm plate	6.57		2		2.75	N/A	
10/29/24	29	21	NONE	KI-CDK9d-32N	5000	2000	3	1x 15-cm plate	6.82		2		2.75	N/A	
10/29/24	29	21	NONE	KB-0742	1500	1500	1	1x 15-cm plate	5.86		2		2.75	N/A	
10/29/24	29	21	NONE	KB-0742	1500	1500	2	1x 15-cm plate	6.27		2		2.75	N/A	
10/29/24	29	21	NONE	KB-0742	1500	1500	3	1x 15-cm plate	6.65		2		2.75	N/A	
10/29/24	29	21	NONE	KB-0742	5000	0	1	1x 15-cm plate							
10/29/24	29	21	NONE	KB-0742	5000	0	2	1x 15-cm plate							
10/29/24	29	21	NONE	KB-0742	5000	0	3	1x 15-cm plate							
10/29/24	29	21	NONE	SY-5609	1000	0	1	1x 15-cm plate							
10/29/24	29	21	NONE	SY-5609	1000	0	2	1x 15-cm plate							
10/29/24	29	21	NONE	SY-5609	1000	0	3	1x 15-cm plate							
10/29/24	29	21	NONE	SY-5609	2500	0	1	1x 15-cm plate							
10/29/24	29	21	NONE	SY-5609	2500	0	2	1x 15-cm plate							
10/29/24	29	21	NONE	SY-5609	2500	0	3	1x 15-cm plate							
10/29/24	29	21	NONE	Senexin B	2000	2000	1	1x 15-cm plate	6.08		2		2.75	N/A	
10/29/24	29	21	NONE	Senexin B	2000	2000	2	1x 15-cm plate	5.99		2		2.75	N/A	
10/29/24	29	21	NONE	Senexin B	2000	2000	3	1x 15-cm plate	5.89		2		2.75	N/A	
10/29/24	29	21	NONE	Senexin B	15000	4000	1	1x 15-cm plate	4.67		2		2.75	N/A	
10/29/24	29	21	NONE	Senexin B	15000	4000	2	1x 15-cm plate	5.72		2		2.75	N/A	
10/29/24	29	21	NONE	Senexin B	15000	4000	3	1x 15-cm plate	5.17		2		2.75	N/A	
10/29/24	29	21	NONE	SEL120	4000	4000	1	1x 15-cm plate	4.26		2		2.75	N/A	
10/29/24	29	21	NONE	SEL120	4000	4000	2	1x 15-cm plate	5.12		2		2.75	N/A	
10/29/24	29	21	NONE	SEL120	4000	4000	3	1x 15-cm plate	4.48		2		2.75	N/A	
10/29/24	29	21	NONE	SEL120	8000	0	1	1x 15-cm plate							No surviving cells
10/29/24	29	21	NONE	SEL120	8000	0	2	1x 15-cm plate							No surviving cells
10/29/24	29	21	NONE	SEL120	8000	0	3	1x 15-cm plate							No surviving cells
10/29/24	29	21	NONE	DMSO	1%	1%	1	1x 15-cm plate	7.86		2		2.75	N/A	
10/29/24	29	21	NONE	DMSO	1%	1%	2	1x 15-cm plate	7.16		2		2.75	N/A	
10/29/24	29	21	NONE	DMSO	1%	1%	3	1x 15-cm plate	8.67		2		2.75	N/A	